

# SHINJect

## **Hub** color

Color coding of the plastic hubs and boxes allows quick and easy identification

### Cannula

The siliconized cannula provides smooth penetration With minimal discomfort

# Lancet point

The ultra-sharp, tri-beveled lancet point minimizes pain and discomfort

### Connection method

Optionally available in Inch thread or metric thread

### **Bevel Indicator**

A bevel indicator on the hub assures easy identification of the lancet bevel position for less insertion pain at the injection site







Office: Room 610-A 6th Floor, Doña Felisa Syjuco

Bldg. Remedios corner Taft Avenue,

Malate, Manila 1004 Philippines Phone: +63 2 8354-2458 / 8526-0896